



Kohli Atraumatic  
**Catheter**



The Urinary Catheter Reinvented  
Better Medicine thru Better Design

# Better Medicine thru Better Design



Conventional Foley Catheter design



Kohli Atraumatic Catheter design

## *When Indwelling Catheters are Part of your Journey*

Introducing a silicone urinary catheter designed to minimize tip contact with the bladder and optimally function for drainage.


The exclusive and patented 5cc pancake balloon tip produces minimal bladder irritation and abrasion<sup>1</sup> while providing improved comfort and drainage.<sup>2</sup>

The Kohli Atraumatic Catheter is designed to meet your long-term catheter needs.

## Ordering Information

Product Code	Description	Unit of Measure
20-100	Kohli Atraumatic Catheter 2-Way, Silicone Urinary Drainage Catheter, 450 mm long, 16 French (5.3 mm diameter)	10/box

The Kohli Atraumatic Catheter is cleared on FDA K172422



 **Kohli Atraumatic  
Catheter**

To place an order, contact Nellie Medical at:

phone 775-800-7300  
fax 844-225-4600  
email [info@kohlicatheter.com](mailto:info@kohlicatheter.com)

[www.KohliCatheter.com](http://www.KohliCatheter.com)

1. J.A. Greenberg, et al, A Preliminary Evaluation of Ovine Bladder Mucosal Damage Associated With 2 Different Indwelling Urinary Catheters, (Journal of Urology, 2017.08.020)

2. L. Hoyte, Bladder Morphology, Emptying, and Residual Volume with 2 Bladder Catheters, (Submission, IUGA 2017)